

# Compound Transfer Spatula

Transfer Powders Easily with this Electronic Spatula!

## Benefits

The handheld electronic spatula collects small quantities of dry compounds for transfer from one vessel to another!

## Features

- Easily transfer very small volumes of powder!
- Best used with standard glass labware



**Integrated with Vial Handling—Coming Soon!**

## Technology

This system replaces the traditional metal laboratory spatula with Transfer Wire. The Transfer Wire, with a controllable electrostatic charge, has the capability to attract, hold, and release powdered samples.

*This compound transfer system utilizes technology developed by INNOVATE Engineering & Design*



**Request a Quote at [biomicrolab.com/contact](http://biomicrolab.com/contact)**